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Fluidsampler instruction.

The Fluid Sampler System will, properly used, give multiple samples from different medias without cross contamination and without the pump being polluted and in need of cleaning.

This is how.

1. Cut a suitable length of tube from the tube roll. Push one end of the tube through the pump head so it barely points out on the other side. Lock the tube to the pump by turning the screw
2. Attach a bottle.
3. Put the other end of the tube into the very place from where you want a sample.
4. Start sucking up a sample by taking just one pump stroke, watching the oil fill the bottle. If needed to get enough quantity, take one or more careful pump stroke.
5. **Very important. Hold the pump with the bottle upright, and don't pump while oil is filling the bottle. In both cases there is a risk for oil to enter the cylinder of the pump, resulting in cross contamination and need for cleaning.**
6. When you have reached desired level, loosen the bottle and the flow will stop.
7. Unscrew the bottle and cap it.
8. Loosen the screw and pull out the tube. Important, if the tube has been emerged in the oil in the bottle, first wipe or cut it off before pulling it out. Otherwise the pump may be polluted.
9. Pull the tube out of the media and dispose of it in an environmental friendly way.
10. You are now ready to take the next sample, without need for cleaning of the pump.

Common errors:

Fierce pumping.

The bottle is not hold upright.

The tube was not clean when it was pulled out.

In above cases there is a risk for polluting of the pump.

If the pump is still being polluted, just unscrew the cylinder, pull back the piston and clean the parts. Be careful not to loose the valve ball or seal ring.