DNA-Marker pUC19/Msp I

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DNA-Marker pUC 19/Msp I

T149.1	50 µg + Probenpuffer
T149.2	4 x 50 µg + Probenpuffer

DNA-Marker pUC19/Msp I

This marker was manufactured by using pUC19-DNA and restriction endonuclease *Msp* I.

Made from natural plasmids: fragment sizes can be validated.

DNA-ladder for analysis of small DNA fragments. The concentration of the individual fragments is equimolar so that rough quantification of the analysed DNA is possible.

The following fragment sizes (in bp) are included:

501, 489, 404, 331, 242, 190, 147, 111, 110, 67, 34, 34, 26.

pUC has a total length of 2686 bp.

The DNA marker is supplied in lyophilised form.

Content:

T149.1 pUC19 /*Msp* I-DNA
0100.1 1 x sample buffer

Application:

Depending on the application, the marker can be absorbed directly in the supplied sterile-filtered 1x sample buffer (ROTI®Load with glycerol) or in TE-buffer (ROTI®Stock 100 x TE, Art. No. 1052.1).

The lyophilized DNA is dissolved for 15 minutes at room temperature in an appropriate volume (see below) of 1x sample buffer ROTI®Load with glycerol (Art. No. 0100.1) under occasional stirring.

1x sample buffer (ROTI®Load with glycerol):

Tris/HCL, pH 7,5	10 mM
Na-acetate	5 mM
EDTA	2 mM
Glycerine	10 % (w/v)
Bromphenol blue	0.02 % (w/v)
Xylene cyanol blue	0.015 % (w/v)

Sample application/concentration:

Standard loading for mini to midi gels per trace is:

- with bands visible in UV-light after Et.-Br.-staining:
0.5 – 1.0 µg in 10 µl = 1 lane
- with detection after Et.-Br.-staining with signal-integrated camera systems:
0.1– 0.5 µg in 10 µl = 1 lane

Storage:

Optimal storage temperature is -20 °C. Please avoid repeated (>10 times) thawing and freezing.

We recommend freezing in aliquots.


Fig.: 1 % agarose (NEEO, Art. No. 2267) + 1 % Synergel™ (Art. No. 0184)

Carl Roth GmbH + Co. KG

Schoemperlenstraße 3-5 • 76185 Karlsruhe
P.O. Box 100121 • 76231 Karlsruhe
Phone: +49 (0) 721/ 5606-0
Fax: +49 (0) 721/ 5606-149
info@carlroth.com • www.carlroth.com ip 07/2021

The company is a limited partnership with headquarters in Karlsruhe, reg. court Mannheim HRA 100055. Roth Chemie GmbH, with headquarters in Karlsruhe, reg. court Mannheim HRB 100428, is the personally liable partner. Managing Director: André Houdelet. Sales tax identification number: DE 143621073.

DNA-Marker pUC 19/Msp I

T149.1 50 µg + sample buffer
T149.2 4 x 50 µg + sample buffer